

3M Infection Prevention Solutions  
Innovation on a Mission

### Steam penetration in steam sterilization processes

Josephus Paulus Clemens Maria van Doornmalen Gomez Hoyos

2011

© 3M 2010. All Rights Reserved.

### Out line

- Steam sterilization conditions
- Loads
- Some physics
- Penetration of steam
- Standards

2

© 3M 2010. All Rights Reserved.

3M Infection Prevention Solutions  
Innovation on a Mission

## Sterilization conditions

3

© 3M 2010. All Rights Reserved.

### Perkins and the Medical Research Council

Perkins 1956<sup>1</sup>

Time [min]	Temp [°C]
2	132
8	125
12	121

<sup>1</sup> Principles and Methods of Sterilization, Perkins JJ, Springfield (IL), 1956  
<sup>2</sup> Working Party on Pressure Steam Sterilizers of the Medical Research Council, Sterilisation by steam under increased pressure, 1959, 273, 425-435

4

© 3M 2010. All Rights Reserved.

### Boiling and egg

5

© 3M 2010. All Rights Reserved.

### Rational Medical Research Council<sup>2</sup>

- Perkins' experiments in water
- Steam sterilization with steam
- Quality of steam not perfect
- Safety margins necessary

Time [min]	Temp [°C]
2 + 1 = 3	132 + 2 = 134
8 + 2 = 10	125 + 2 = 127
12 + 3 = 15	121 + 0 = 121


<sup>2</sup> Working Party on Pressure Steam Sterilizers of the Medical Research Council, Sterilisation by steam under increased pressure, 1959, 273, 425-435

6

© 3M 2010. All Rights Reserved.

### Steam sterilization conditions


- On all surfaces 100 % or saturated steam
- Predetermined temperature
  - E.g. 134 °C and 121 °C
- Predetermined time
  - E.g. 134 °C for 3 minutes and 121 °C for 15 minutes



7

3M Infection Prevention Solutions  
Innovation with a Mission


# Loads




8

### Surfaces in porous loads

- Textiles
- Bandages
- Cotton-wool





All have 'inner' surfaces



9

### Instrument surfaces


- Surfaces of instruments

10

### Change of loads in hospitals

- Porous loads are replaced by disposables or single use
  - E.g., cotton drapes replaced by non woven and crepe
- Instruments become more complex
  - E.g., More Minimal Invasive Surgery result in more lumened instruments



11

### Capillary suction in porous loads



Decreasing diameter capillary suction



12

3M Infection Prevention Solutions  
Innovation  
with a Mission

# Sterilization

**3M**

13 © 3M 2010. All Rights Reserved.

Sterilization is a combination of

- Sterilizer
- Process
- Load
- Loading pattern
- Wrapping (Micro Biological Barrier)


Changes in one may have changes on the end result

**3M**


14 © 3M 2010. All Rights Reserved.

Example

Change:  
Load from non hollow instruments to hollow instruments




Action:  
To establish steam sterilization conditions on all surfaces  
Process from gravity to fractionated vacuum process



**3M**

15 © 3M 2010. All Rights Reserved.

Note



Because of capillary working of porous loads:  
It is possible to sterilize porous loads in super atmospheric processes

Disadvantages

- Less reproducible
- Longer process time

**3M**


16 © 3M 2010. All Rights Reserved.

Steam Sterilization processes

**Sterilization condition established by:**  
Sterilization condition Must be met at the end of the defined sterilization phase: temperature and time

Phase 1: Conditioning  
Phase 2: Sterilization  
Phase 3: Drying and pressure equilibration

Convection  
Diffusion  
Condensation



All kind of different air replacement methods  
No standards method

**3M**

17 © 3M 2010. All Rights Reserved.

3M Infection Prevention Solutions  
Innovation  
with a Mission

# Convection

**3M**

18 © 3M 2010. All Rights Reserved.

**We control the flow**

Convection

Ceiling

Room

Heater

3M

19

Only convection

Process

Pressure

Time

Steam inlet

Pumping of gas mix

Sterilizer chamber

3M

20

Only convection and tube in sterilizer chamber

Sterilizer chamber

Process

Compressing and decompressing the gas

Convection can be imposed

Time

Steam inlet

Pumping of gas mix

Interface between air and steam  
Act like a piston in a cylinder

3M

21

3M Infection Prevention Solutions  
Innovation  
with a Mission

Diffusion

3M

22

Diffusion

**Diffusion needs time**

Diffusion cannot be imposed

Homogenous distribution of the gas

Bottle of Perfume

Room

3M

23

Diffusion in lumen

Air

Steam

Air diffuses into the steam  
and  
Steam diffuses into the air

Diffusion cannot be imposed.  
It is dictated by nature

3M

24

### Diffusion influences the air removal time

Fast pulsing

Slow pulsing

Process

Pressure

Time

More time for diffusion

In steam sterilization processes diffusion is slower than convection

With fast pulsing **NO** time for diffusion

**3M**

3M Infection Prevention Solutions  
Innovation  
with a Mission

# Condensation

**3M**

### Condensation

- Basic of condensation are understood
  - Basically:
    - Steam condenses on colder surfaces
    - Latent energy transferred to surface

Condensation steam  
Volume reduction in the order  
of 1800 times

Cold surface

**3M**

### Example steam needed to warm up stainless steel

- 1 kg stainless steel instruments
  - Warm up for 24 to 134 °C (difference = 110 °C)
  - Heat capacity stainless steel  $\approx 4.600 \text{ J} / (\text{kg} \cdot ^\circ\text{C})$
  - Needs  $1 \text{ kg} \times 4.600 \text{ J} / (\text{kg} \cdot ^\circ\text{C}) \times 110 \text{ }^\circ\text{C} = 50.600 \text{ J} = 50.6 \text{ kJ}$
- Steam at 134 °C
  - Energy of steam  $\approx 2700 \text{ kJ/kg}$
- Steam needed to warm up instruments:
  - $(50.6 \text{ kJ}) / (2700 \text{ kJ/kg}) = 0.019 \text{ kg}$
  - 0.019 kg steam equals about  $(0.019 \text{ kg}) / (0.60 \text{ m}^3/\text{kg}) = 0.032 \text{ m}^3$
  - 0.032 m<sup>3</sup> = 32 l steam

**3M**

### Back to lumen

Temperature lumen equals temperature steam

Temperature lumen lower temperature steam

Steam condenses immediately on wall

Volume reduction of about 1700 times

More steam supplied to equalize pressure

Process continues until wall is warmer or equal to steam temperature

**3M**




### Physical difference steam and Air (NCGs)

Condense forming

Steam keeps streaming in  
Energy is transported in

Air shrinks a little bit, no flow  
Energy transport is slow

**3M**



### Conductivity and velocity



$$\frac{\text{Distance difference}}{\text{Time difference}} \text{ in } \frac{\text{km}}{\text{hour}}$$




$$\frac{\text{Temperature difference}}{\text{Time difference}} = \frac{\Delta T}{\Delta t} = \frac{dT}{dt} \text{ in } \frac{^\circ\text{C}}{\text{s}}$$

Conductivity =  $\frac{dT}{dt}$  X with some properties of the material

in  $\frac{\text{W}}{\text{m}^\circ\text{C}}$

34 © 3M 2010. All Rights Reserved.





- ### Causes for presents of NCGs resulting fails
- Insufficient air removal, e.g.,
    - Not deep enough vacuum
    - Not high enough steam injections
    - Too fast pulsing
  - Leak in vessel, pipe work, valves or gaskets
  - NCGs in steam
- 3M © 2010. All Rights Reserved.

- ### Consequently
- By measuring NCGs:
    - Steam sterilization conditions can be checked
    - State of the art technology
- 3M © 2010. All Rights Reserved.

