


The Removal and Inactivation of Prions



December 6th to 8th, 2007
Bydgoszcz, Poland
Peter Newson
Marketing Manager -EMEA

STERIS




■ Agenda


1. *Introduction*
 1. *Subject*
 2. *Risks*
 3. *Standards*

1. *Prion introduction*

3. *Cleaning*
 1. *Decontamination*
 2. *Prion inactivation*



■ Target



Inactivation of prions on infectious proteins

■ The subject



- **Contaminated surgical instruments need to be cleaned to remove soil and infectious items and proteins**

■ Soil-A Dirty Subject



- **Implications of inefficient cleaning**
 - Visually unacceptable devices
 - Ineffective disinfection and sterilization
 - Device damage
 - Adverse patient outcomes
 - Accidental prion transmission

■ Increased Attention

- **Visually unacceptable devices**
 - Cancelled procedures
- **Device damage**
 - Dried on soil / salt deposits
- **End toxins**
- **TASS**
 - Toxic Anterior Segment Syndrome
 - Outbreak in the USA
 - Cataract surgery
 - Inadequate cleaning and/or water residuals



■ What has Changed since 2000 ?

- **The Science of cleaning**
- **Understanding of “Prions”**
- **Revolution in Prion Decontamination**
 - New Investigations
 - New products

- **New Regulatory Focus**

■ What has changed.....

.....and reasons for the development of new and more effective decontamination methods:

- **Patient risks**
- **New Standards**
- **Increased attention**

■ Patient risks

- **Inadequate disinfection and sterilization**
- **Toxin risks**
- **Prion diseases**
 - Other protein disease concerns



■ Standards for decontamination

- EN/ISO 15883 series
 - Published
 - *Part 1: General requirements, definitions and tests*
 - *Part 2: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, hollowware, utensils, glassware, etc.*
 - *Part 3: Requirements and tests for washer-disinfectors employing thermal disinfection for human waste containers*
 - In final draft (FDIS)
 - *Part 4: Requirements and tests for washer-disinfectors employing chemical disinfection for thermo-labile endoscopes*
 - Technical Specification (TS)
 - *Part 5: Test soils and methods for demonstrating cleaning efficacy of washer-disinfectors*




■ Standards for sterilization

- EN/ISO 17665, *Sterilization of medical devices- Information to be provided by the manufacturer for the processing of resterilizable medical devices*
- Others
 - Biocides directive (98/8/EC; April, 1998)
 - Detergents Directive (648/2004; October, 2005)
 - REACH (Com(03) 644)
 - ‘Registration Evaluation and Authorisation of Chemicals’




■ Prion introduction

- What are 'Prions'?
 - Identified as the causative agents for a group of central nervous system diseases
 - TSEs (Transmissible Spongiform Encephalopathies)
 - vCJD, CJD (Creutzfeldt-Jacob disease)
 - *Still debated!*
 - Proteins
 - Appear to be devoid of nucleic acid



■ Prion introduction


Prions and protein folding: The cause of vCJD appears to be the presence of an abnormal forms of normal protein, known as a **prion**. This structure makes them hard to be removed from the instruments.



- Accumulate in the body and Causing cell death (brain).
- No more cell growth


■ Prion decontamination

- Intrinsic resistance
 - Prions demonstrate resistance to routine methods of decontamination and sterilization
 - Prions are proteins, not microorganisms




■ Accidental transmission in practice

Source of CJD




Case 1



Case 2



Assay



Implantation



Patient (female)
Age: 69

Exposure Sept. 1974

Patient (female)
Age: 23

Exposure Nov. 1974
Death 20 month later

Patient (male)
Age: 17

Exposure Dec. 1974
Death 16 month later

Chimpanzee
Age: 17

Exposure Dec. 1977
Death 20 month later



1-2 new diseases per 1 Mio. People a Year / Global

■ STERIS Research Facility

- Fontenay-aux-Roses, Paris
 - Laboratory facilities
 - Microbiology (Biosafety 1-3), histology and biochemical laboratories
 - Atypical pathogen research
 - Including Prions



In Vivo and ***In vitro*** sets where carried out
In collaboration with STERIS

The Removal and Inactivation of Prions

CLEANING

STERIS



■ The science around cleaning

- The chemistry of soil
- How much and what is present
- Detection of soil
- Cleaning indicators
- ‘New’ cleaning methods
-How clean is clean?



■ **Concerns with cleaning**

- Infectivity is difficult to remove from instrument surfaces
- Drying of soil increases risk
- Improper handling increases risk
- Some standard cleaning processes may also increase risk
- Impractical recommendations

■ Considerations

- Efficacy
 - Data is conflicting
 - Testing on surfaces
- Compatibility
 - Surface damage
- Safety
 - Handling



WHO process: Emerge instruments in a High concentration of Sodium Hydroxide

■ The cleaning process


Where should reprocessing begin?

- Point of use (OR)
- WHO, NICE, AORN Guidelines-
recommend instruments should not be
allowed to dry

Instrument transportation gel / pre-soak

■ Important Detergent Concepts

- **Regulations**
 - Claims
- **Safety**
- **Process**
 - Contact Method, Time, Temperature, Concentration
 - Water quality
- **Formulation**
 - A combination of ingredients, including active and inert ingredients, into a product for its intended use



■ New Cleaning Technology

- Enzymatic formulations
 - Genencor
 - Alkaline proteases
 - Effect of formulation
- Advances in alkaline cleaning
- Physico-chemical cleaning
 - e.g., plasma

■ Cleaning performance

- Cleaning performance test on static emulsion

<i>Enzymatic Methods</i>	<i>„Log“ Reduction</i>
Klenzyme	~4.5
Enzyme Cleaner 2	~1
Porous load autoclaving 15 min at 121°C	~5.5
Klenzyme + autoclaving	~6.5
Enzyme Cleaner 2 + autoclaving	~3.0

<i>Alkaline Methods</i>	<i>„Log“ Reduction</i>
Hamo50 (1%, 65°C; 10 min.)	~3
Hamo50 (1%, 70°C; 10 min.)	~6.5
CIP150 (1,6%, 43°C; 15 min.)	>7
Hamo100 (1%, 65°C; 10 min.)	>7
(0,8%; 43°C; 7,7 min.)	>7
(0,8%; 43°C; 7,7 min.); Plastic wires	>7

Not all enzymatic & alkaline cleaners perform in the same way



Detergent formulations are critical for prion inactivation & removal

■ Two in one Detergent

- Outstanding Cleaner
- Removes and Inactivates Prions

■ Key product components

- Potassium Hydroxide (KOH)
- Surfactants
- Chelating Agents

■ Key product components

Why Potassium Hydroxide (KOH) over Sodium Hydroxide (NaOH)?

- Both increase the solubility of soils
- Potassium ions and salts are more soluble in water than Sodium ions and salts
- No need to neutralize with acid
- **Better Rinsing Ability!**

■ Key product components

- **Surfactants- Surface Active Agent**
 - Penetrate soil and surface irregularities
 - Ability to displace particles



No Surfactants



Surfactant A



Surfactant B

■ Key product components

Chelating Agents

- Sequester hard water ions
- Keeps them in a suspension
- Prevents them re-depositing on surface

■ Prion Decontamination Cycle

• Manual Soak with HAMO100 PID

- Fill basin with water to a temp of 43°C (keep warm)
- Add **0.8%** Hamo 100 PID
- Submerge instruments in water (5cm above instruments)
- Typical soak time = **7.5 min**
- Rinse and then pass through normal process



■ Prion Decontamination Cycle

- **HAMO 100 PID usage in washer/disinfector**


Phase	Temp	Time	Concentration	Comment
Pre- Wash	$\leq 43\text{ }^{\circ}\text{C}$	1 – 2.5 min	0.2%	Dosed pre-wash can be used when instruments have been allowed to dry
Wash	$\geq 43\text{ }^{\circ}\text{C}$	4 – 7.5 min	0.2 – 0.8%	Concentration depends on water quality

■ How does the formulation act

- **Removes** prions & bioburden from surfaces
- Keeps them **suspended**
- **Inactivates** through **hydrolysis** & allows suspended residues to be **freely rinsed** without re-adhering to the surface
- **Substrate Compatibility** (Instrument materials)

■ Hydrolysis as inactivating mechanism



- Breaks down the peptide bonds through the introduction of water between the atoms of the bond, thus inactivating the infectivity of the prion



■ Degradation with Hamo 100

- >6 log reduction in prion titers against multiple strains of prions including CJD & vCJD
- Total cleaning of ISOEN15883 part 5 test soils incl. Edinburgh & blood soils

■ Degradation with Hamo 100



■ Conclusion

- Prions can be transmitted on device surfaces
- Methods have been developed to verify the effectiveness of decontamination methods
- Cleaning alone can be effective
 - Choice of cleaner and process is important
- Disinfection and sterilisation can be effective
 - Choice of biocide and process is important
 - New methods of prion decontamination have been identified

The Removal and Inactivation of Prions



Peter Newson

