

What is the place for Hydrogen peroxide plasma sterilisation in hospital sterile supply?

Walter KOLLER, Univ.Prof., MD

Head, Clinical Division for Hospital Hygiene
Clin.Institute for Hygiene and Med. Microbiology

Medical University Vienna

Background



- Growing interest in easy-to-operate sterilisers for complex medical devices (MD)
- Simplicity and reliability claimed for Sterrad™ hydrogen peroxide vapour (HPV) sterilisers
- Conflicting results from different investigator teams
- German/Austrian multicenter study warranted restrictions for HPV-steriliser use in clinical settings

Study protocol (1): authors

- Three STERRAD 100S sterilisers tested in parallel with the same protocol in Vienna/A (two centres ^{1, 2}) and Schwerin/Ge ³

¹ KOLLER W, Clin.Div. Hospital Hygiene, Medical University Vienna, Austria

² GETREUER H, Dept. for Hospital affairs, Vienna County, Austria

³ WERNER H P, Hyg-Cen, Schwerin, and KRAMER A, Dept. Hygiene and Environmental Medicine, Ernst-Moritz-Arndt University, Greifswald, Germany

Study protocol (2): micro lab

- 10^7 /ml suspension of *G. stearothermophilus* dried onto carriers and exposed to HPV sterilisation cycle
- Spore recovery from carriers by vortexing in CSB and plating in serial dilutions on CSA
- Spore killing effect expressed as log reduction calculated from exposed and unexposed carriers
- Each variation tested in triplicate

Study protocol (3): variables



- Spore suspension media
- Carrier types
- Carrier wrapping
- Steriliser load size
- Steriliser cycle length
- Carrier position in steriliser chamber

Legends to figures ⁽¹⁾

LOG-REDuctions as
calculated after HPV
exposure in Sterrad

(LOG-RED = $\log \text{cfu}_{\text{before}} - \log \text{cfu}_{\text{after}}$)

SAL 10^{-6}

= sterility assurance level

LABoratories/
Investigators (city):

- **koller** (Vienna)
- **werner/Kramer**
(Schwerin/Greifswald)
- **getreuer** (Vienna)

Legends to figures ⁽²⁾

SPORe SUSPensions:

G.stearothermophilus-spores ($\sim 10^6$ /carrier) in:

napep isotonic saline + peptone

adest dest. water

Legends to figures ⁽³⁾



EXPOSITION of spores: carriers/devices
(wrap- and *carrier-type*):


PEhelix polyethylene-helix (*fine mesh metal*)

Methelix metal-helix (*fine mesh metal*)

Tyvekfol Tyvek-wrap (*fine mesh metal*)



Wraps for fine mesh metal spore-carriers

PEhelix




Tyvekfol

Methelix







Var's: **spore suspensions** and **carriers/expositions**




General effect of position in dose/response curves on variances




Effect of **steriliser load size** (spore suspension **adest**)



Effect of **cycle length** (spore suspension **adest**)



Effect of **spore carrier position** in steriliser (spore susp **adest**)



Conclusions from study 2000

Sterrad™ 100S only acceptable with significant limitations:

- Absence of any organic or crystalline mini-residues on surfaces
- Absence of long fine bore channels esp. in metallic structures
- Steriliser load below critical limit
- Absence of cotton/cellulose and other organic components in MD and wrap

Control and documentation of critical process parameters urged!

2007 status

New generation of STERRAD 100-S, NX and 200 systems equipped with process control and documentation tools which allow validation.

In October 2006, a **limited approval was issued by the Austrian ministry of health (AMH)** for these STERRAD products. **AMH Limitations** for Sterrad™ use in human medicine comprise:

- Compliance with EN ISO 14937:2000 mandatory for **use in medical settings** (notably documentation and control of all relevant process parameters with suitable and calibrated instruments)
- optimised and validated precleaning mandatory

Limitations (cont.)

- exclusion of any MD which can be autoclaved
- exclusion of (invasive) MD equipped with hardly accessible fine bores, joints, seams, grooves, lamellae or grills
- STERRAD not accepted as alternative to EO- or LTF-sterilisation unless suitability of and compatibility with MD in the given clinical setting proven and on site validation successfully performed
- Specialised instrument wraps (Tyvec) and sealing equipment as recommended by mfct. to be used in any case.

Recommendations

- New STERRAD 100-S, NX and 200 systems recommended for MD which
 - require sterilisation and
 - are not intended for invasive use and
 - do not require special manipulations in the process steps preceding sterilisation

= exclusively for highly specialised procedures and equipment !
- Specialised instrument wraps (Tyvec) and sealing equipment must be budgeted when comparing with other LT sterilising systems

Conclusions

- HPV can act sporocidal at moderate temperatures (40 - 60 °C) when high level cleanness and uninhibited access of HPV is guaranteed.
- New STERRAD 100-S, NX and 200 systems are equipped with process control and documentation tools which allow validated sterilising processes,
- approved by the Austrian ministry of health as niche products for sterilisation of highly selected MDs under strictly controlled conditions (EN ISO 14937:2000).
- Specialised instrument wraps (Tyvec) and sealing equipment are required.